

# TECHNICAL SPECIFICATIONS • C300 PS

## MEASUREMENT AND PERFORMANCE

Length	.45 in
Minimum transport length	.34.5 in
Width	.24.25 - 33 in
Height (floor to top of standard back)	.41 in
Seat to floor height	.17 in.
Seat width	.16/18/20/22 in
Seat depth	.16-22 in
Backrest height	.18.5/22 in
Armrest height	.8.5-13.5 in
Armrest width	.2/3/4 in
Backrest angle - adjustable	.85°-115°
Leg rest angle	.70°
Seat elevator	.8 in
Max. user weight	.300 lbs
Weight incl. batteries	.285 lbs
Range	.16 miles
Max. speed (forward)	.5 mph
Turning space   80°	.40 in
Ground clearance	.3 in
Vertical obstacle climb	.2.5 in
Uphill angle	.10°
Downhill angle	.6°
Front tire size	.300 x 8
Rear tire size	.2.10 x 65
Battery capacity	.2x60AH
Battery needed (x2)	.GR 34
Charging time	.approx. 8 hours
Electronics available	.VSI

- Super-low seat to floor height of 17 inches
- Increased ground clearance to 3 inches
- Front and rear battery access
- Adjustable suspension
- Virtually limitless angle adjustments for footplate
- Fixed or adjustable height armrests, both flip-back and are removable
- Universal 1 inch tubing to support aftermarket options and accessories

## STANDARD OPTIONS

Color options: Vintage Wine, Crystal Blue  
Choice of seat width/depth: 16/18/20/22x16-22 in  
Choice of fixed or adjustable armrest height

## STANDARD ACCESSORIES

Optional joystick knobs  
Positioning belt  
Back Pack Hooks  
Adjustable head support



**Permobil, Inc.**  
6961 Eastgate Blvd.  
Lebanon, TN 37090  
Tel: 800-736-0925 Fax: 800-231-3256  
[www.permobil.com](http://www.permobil.com)