

Aluminum Commode is Lightweight and Durable

NEW!



MDS89664A

- High-strength, aluminum frame prevents rusting.
- Durable, plastic seat with lid.
- Plastic armrests for additional, patient comfort.
- Adjustable legs allow the seat to be raised or lowered in 1" increments to sit from 16" to 20".
- Includes 10-quart capacity, commode bucket with cover and splash shield.
- Non-skid, non-marring, rubber tips provide extra safety.
- 300 lb. weight capacity.

Item #	Description	Base		Width		Seat		Seat		Height		Weight	Pkg.	HCPCS
		Width	Depth	Between	Arms	Depth	Width	Overall	Min./Max.	Min./Max.	Capacity			
MDS89664A	Aluminum	25½"	16½"	18½"	14"	15"	16"-20"	23"-27"	300 lbs.	4/cs	E0163			

Deluxe, Steel Commodes Easily Convert for Use as a Commode, Raised Toilet Seat or a Toilet Safety Frame




MDS89664



MDS89664XW

- Back rail removes quickly and easily.
- Features a durable, plastic seat with lid.
- Plastic armrests for additional patient comfort.
- Adjustable legs allow seat to be raised or lowered in 1" increments to sit from 16" to 22½" standard, 19½" to 23½" bariatric.
- Includes 10-quart capacity, commode bucket w/cover & splash shield.
- Non-skid, non-marring, rubber tips provide extra safety.



Item #	Description	Base		Width		Seat		Seat		Height		Weight	Pkg.	HCPCS
		Width	Depth	Between	Arms	Depth	Width	Overall	Min./Max.	Min./Max.	Capacity			
MDS89664	Enamel, 4-in-1	24½"	18"	19"	14"	15"	16"-22½"	23½"-30"	300 lbs.	4/cs	E0163			
MDS89664XW	Bariatric Commode	28"	20½"	24"	14"	15"	19½"-23½"	29½"-32½"	650 lbs.	each	E0168			

Drop-Arm Steel Commodes Makes Patient Transfers Easy



MDS89668



MDS89668XW

- Angled front is comfortable for cast users.
- Single-handed, arm release.
- May be used over standard toilet.
- Lowered arm supports weight during transfer.
- Chrome-plated, steel frame.
- Five-way pail removal.
- 12-quart pail w/cover included.



Item #	Description	Base		Width		Seat		Seat		Height		Weight	Pkg.	HCPCS
		Width	Depth	Between	Arms	Depth	Width	Overall	Min./Max.	Min./Max.	Capacity			
MDS89668	Drop Arm Steel	20¼"	20½"	17½"	14"	15"	17"-21"	28¾"-33¾"	250 lbs.	each	E0165			
MDS89668XW	Bariatric Commode	25"	19"	24"	23"	18"	19¼"-23¼"	25½"-30½"	850 lbs.	each	E0168			

Universal Fit Commode Bucket

- Constructed of high-impact, white polyethylene.
- Universal - Fits all Medline commodes.
- Complete with handle and lid.
- Diameter of bucket is 11".
- 12-quart/ 2-gallon capacity.



MDS80306B

Commode Accessories

Item #	Description	Pkg.
MDS89664SL	Seat & lid for MDS89664	4 sets/cs
MDS89664SLXW	Seat & lid for MDS89664XW	1 set
MDS89664T	Set of four tips for MDS89664 & MDS89664XW	4 tips/set
MDS89664SGC	Set of four splash guards	4 splash guards/set

Item #	Description	Pkg.	HCPCS
MDS80306B	Commode Bucket	6/cs	E0167