

Locking, Raised Toilet Seats Help Prevent Shifting While Sitting or Transferring for Optimum Patient Safety and Comfort



- Width adjustable, 2 position arms eliminate the need for toilet safety rails, making it easier for the patient.
- Front clamping mechanism ensures positive locking onto toilet. Helps prevent shifting of toilet seat while sitting or transferring.
- Fits most toilet bowls from 11½"-15" long.
- Wide, contoured surface provides maximum user comfort.
- Padded support arms offer user a secure hand support and maximum comfort when using toilet seat.
- Raises seat level by 5½".
- 300 lb. weight capacity.

Item #	Description	Height	Seat	Seat	Wt. Cpcty.	Pkg.	HCPCS
		Min./Max.	Width	Depth			
MDS80314	Raised Toilet Seat w/o Arms	5½"	16⅞"	17⅞"	300 lbs.	3/cs	E0244
MDS80316	Raised Toilet Seat with Arms	5½"	16⅞"	17⅞"	300 lbs.	3/cs	E0244

Universal Raised Toilet Seat Fits Securely Over Any Toilet



- Seat fits securely over any toilet bowl.
- Vinyl-covered brackets allow height adjustments from 4¾" to 6¾" above standard seat and prevent chipping or scratching of porcelain.
- Almond plastic seat and plastic splash guard make cleaning easy.
- 300 lb. weight capacity.



MDS80315

Item #	Description	Height	Seat	Seat	Wt. Cpcty.	Pkg.	HCPCS
		Min./Max.	Width	Depth			
MDS80315	Universal Raised Toilet Seat	4¾"-6¾"	15"	15½"	300 lbs.	each	E0244



MDS80318RW

Molded Plastic Raised Toilet Seat is Contoured for Patient Comfort



- Contoured for patient comfort.
- Heavy-duty, white polyethylene seat provides excellent support.
- Raises seat level by 4½".
- 300 lb. weight capacity.

Item #	Description	Height	Seat	Seat	Wt. Cpcty.	Pkg.	HCPCS
		Min./Max.	Width	Depth			
MDS80318RW	Plastic Raised Toilet Seat	4½"	15"	16"	300 lbs.	3/cs	E0244