

EXCEL 2000

SPECIFICATIONS

Weight Capacity	300 lbs
Seat Width	18", 20"
Seat Depth	16"
Standard Chart Pocket	yes
Footplates	composite
Dual Axle	yes
HCPCS Code	K0001, K0002
For full specifications see pages 28-29	

ACCESSORIES:

O ₂ Holder	MDS85181U
I.V. Pole	MDS85183
O ₂ / I.V. Combo	MDS85190
Anti-Fold Device	MDS85196
Overhead Anti-Theft Device	MDS85197A
Anti-Tip Device	MDS85189



MDS806300D

FEATURES:

- Carbon steel frame with chip-resistant chrome plating.
- Powder-coated cross bars and clothing guard.
- Tough, embossed upholstery.
- Threaded seat screw inserts reduce stripping from upholstery changes.
- Mag-style composite wheels are extremely durable and easy to clean.
- Heavy-duty fork and welded fork stems minimize wobbling.
- Sealed bearings in fork housing, rear wheels and front wheels keep dirt out and are maintenance free, adding durability.
- Padded armrests.
- Swing-away elevating legrest models feature padded calf pads.
- Elevating legrest models feature a notched, stainless steel ratchet bar to lock the legrests securely in place.

Item #	Description	Available Upholstery Color(s)
MDS806100*	18" Permanent full-length arms, fixed footrests	Ruby, Navy
MDS806150*	18" Permanent full-length arms, swing-away detachable footrests	Ruby, Navy
MDS806200*	18" Permanent full-length arms, swing-away detachable elevating legrests	Ruby, Navy
MDS806250*	18" Removable desk-length arms, swing-away detachable footrests	Ruby, Navy
MDS806300*	18" Removable desk-length arms, swing-away detachable elevating legrests	Ruby, Navy
MDS806400	20" Removable desk-length arms, swing-away detachable footrests	Black
MDS806450	20" Removable desk-length arms, swing-away detachable elevating legrests	Black
MDS806100D	18" Permanent full-length arms, fixed footrests	Black
MDS806150D	18" Permanent full-length arms, swing-away detachable footrests	Black
MDS806200D	18" Permanent full-length arms, swing-away detachable elevating legrests	Black
MDS806250D	18" Removable desk-length arms, swing-away detachable footrests	Black
MDS806250DFLA	18" Removable full-length arms, swing-away detachable footrests	Black
MDS806300D	18" Removable desk-length arms, swing-away detachable elevating legrests	Black
MDS806300DFLA	18" Removable full-length arms, swing-away detachable elevating legrests	Black

COLOR OPTIONS



*To order, add the three letter color code designating color to the end of the item number.
Ruby—RBY; Navy—NVY All other chairs listed are standard black upholstery