

EXCEL K1 Basic

SPECIFICATIONS

Weight Capacity	250 lbs
Seat Width	18"
Seat Depth	16"
Standard Chart Pocket	no
Footplates	composite
Dual Axle	yes
HCPCS Code	pending
For full specifications see pages	28-29

ACCESSORIES :

O ₂ Holder	MDS85181U
I.V. Pole	MDS85183
O ₂ / I.V. Combo	MDS85190
Anti-Fold Device	MDS85196
Overhead Anti-Theft Device	MDS85197A
Anti-Tip Device	MDS85189



MDS806300EE

FEATURES :

- Durable, tig-welded gray frame.
- Hemi-height adjustable seat and front casters.
- 250 lb. weight capacity.
- Comfortable nylon upholstery.
- Smooth-rolling latex tires.
- Lifetime frame warranty; 1 year on all other parts.
- K1, K2 coding (pending).

Item #	Description
MDS806250NEE	16" Desk-length arms, swing-away detachable footrests
MDS806150EE	18" Permanent full-length arms, swing-away detachable footrests
MDS806250EE	18" Desk-length arms, swing-away detachable footrests
MDS806300EE	18" Desk-length arms, swing-away detachable elevating legrests
Upholstery Color: Black	