

# EXCEL Kidz

## SPECIFICATIONS

Weight Capacity	250 lbs
Seat Width	12", 14"
Seat Depth	12"
Standard Chart Pocket	no
Footplates	composite
Dual Axle	yes
HCPCS Code	pending
For full specifications see pages	28-29

## ACCESSORIES :

O <sub>2</sub> Holder	MDS85181U
I.V. Pole	MDS85183
O <sub>2</sub> / I.V. Combo	MDS85190
Anti-Fold Device	MDS85196K
Overhead Anti-Theft Device	MDS85197K



MDS806120PEDE

## FEATURES :

- Tough nylon upholstery
- Swing-back arms.
- Cast aluminum seat cradles (instead of plastic) for durability.
- Standard dual axle allows for easy adjustment of seat height to hemi level.
- Treaded tires for safe riding.
- Padded and upholstered armrests are standard features.
- Anti-Tippers included.
- Smiley face on seat back.
- Telescoping handles for caregiver **(COMING SOON!)**

Item #	Description	Frame Color
<b>MDS806120PED</b>	12" Swing-back desk-length arms, swing-away footrests	Yellow
<b>MDS806120PEDE</b>	12" Swing-back desk-length arms, elevating footrests	Yellow
<b>MDS806140PED</b>	14" Swing-back desk-length arms, swing-away footrests	Yellow
<b>MDS806140PEDE</b>	14" Swing-back desk-length arms, elevating footrests	Yellow

Upholstery Color: Black with Yellow Smiley Face on Seat Back