

# EXCEL Superlight Freedom Transport Chair

## SPECIFICATIONS

Weight Capacity	300 lbs
Seat Width	19"
Seat Depth	16"
Standard Chart Pocket	no
Footplates	composite
Dual Axle	no
HCPC Code	pending
For full specifications see pages	28-29

## ACCESSORIES :

O <sub>2</sub> Holder	MDS85181FT*
I.V. Pole	MDS85183FT*
O <sub>2</sub> / I.V. Combo	MDS85190FT*
Anti-Fold Anti-Theft Device	MDS85196
*Anti-Tip Device	MDS85189FT
Recommended for use with Anti-Tip Device	



## FEATURES :

- Weighs less than 15 pounds (25% less than competitors).
- Nylon upholstery.
- Cup holder standard.
- Backside features bottle holder, key ring, cell phone pocket, identification holder, and pocket for billfold or iPod.
- Back folds down for easy storage and transport.
- Convenient carrying handles on seat.



Backside of chair



Optional carrying case with shoulder strap (MDSCHAIRCASE)

Item #	Description	Frame Color
<b>MDS808200SLR</b>	Permanent desk-length arms, swing-away detachable footrests	Red
<b>MDS808200SLB</b>	Permanent desk-length arms, swing-away detachable footrests	Blue
Upholstery Color: Black		