

TECHNICAL SPECIFICATIONS • C500 CORPUS

MEASUREMENT AND PERFORMANCE

Length	.38 in
Minimum transport length	.33 in
Width	.25.75 - 28.5 in
Height (floor to top of standard back)	.46 in
Seat to floor height	.19.75 in
Seat width	.17/19/21/23 in
Seat depth	.14-22 in
Backrest height	.21/25.5 in
Armrest height	.5-12.5 in
Armrest width	.3-4 in
Backrest angle	.90°-175°
Leg rest angle	.90°-170°
Seat elevator	.8 in
Max. user weight (310 lbs. with HD package)	.250 lbs
Weight incl. batteries	.350 lbs
Range	.19-25 miles
Max. speed (forward)	.5 mph
Turning space 180°	.44 in
Ground clearance	.2.75 in
Vertical obstacle climb	.3 in
Uphill angle	.12°
Downhill angle	.12°
Front tire size	.300 x 8
Rear tire size	.2.10 x 65
Motor type	.4 pole
Battery capacity	.2x60AH
Battery needed (x2)	.GR 24
Charging time	.approx. 8 hours
Electronics available	.Pilot+

- Speed package available with a top speed of 7.5 mph
- Ability to maneuver in tight living spaces with a 26" turning radius
- Excellent option for indoor and outdoor users
- Four pole motors providing power and torque to power through rough terrain
- Adjustable independent suspension providing a smooth ride
- Ergonomic back with built-in, adjustable lumbar and trunk supports
- Availability of additional seat components, such as lateral and thigh supports
- Power functions may include: tilt, recline, seat elevator, legrest elevation and many custom powered accessories

STANDARD OPTIONS

Color options: Vintage Wine, Crystal Blue

Choice of seat width/depth: 17x16/18/20/22 in, 19x18/20/22 in, 21x18/20/22 in, 23x18/20/22

Choice of backrest width/height: 14x21/25.5 in, 16x21/25.5/29.5 in, 18x21/25.5 in

Choice of armrest pad width/length: 4x10/12/14/16 in

STANDARD ACCESSORIES

Optional joystick knobs

Swing away trunk supports

Lateral Wedges

Thigh Supports

Push Handles

Positioning belt

Removable leatherette headrest

Back Pack Hooks



Permobil, Inc.
6961 Eastgate Blvd.
Lebanon, TN 37090
Tel: 800-736-0925 Fax: 800-231-3256
www.permobil.com