

MEASUREMENTS AND PERFORMANCE

Length	38"
Width	22"
Height	28"/30"/34"
Transportable size	(L/W/H) 29" / 22" / 26"
Weight including batteries	190 lbs.
Maximum user weight	110 lbs.

BASE

Base shroud color	yellow/red
Battery type	M40 gel cell (recommended - M40 12 SLD G) ¹
Charging time	6-8 hours (approx.)
Driving range	16 miles ²
Maximum speed	4.5 mph
Turning radius	22"
Turning space, 180°	35"
Maximum curb height	2.5"

SEAT

Choice of upholstery	leatherette, fabric
Seat height	(lowered – elevated) 15" – 23"
Seat depth	(S/M/L) 10" / 12" / 14"
Seat width	(S/M/L) 10" / 12" / 14"
Back rest height	(S/M/L) 14" / 16" / 20"
Back rest width	(S/M/L) 10" / 12" / 14"
Back rest angle	(manual recline) 90° – 165°
Seat tilt angle	.0° - 25° / -5° - 20°
Arm rest height	5" – 10"
Foot rest distance from seat pan	5" – 13"
Leg support angle	(manual adj.) 90° – 140°
Foot plate in two sizes	.S and L



- Pediatric wheelchair for toddlers to young teens
- Unique low seat to floor height combined with the small footprint allows for seamless interaction with peers.
- Equipped with flexible Mini-Flex seat, which lets the child grow in the seat. Change the seat plates and cushions when you need a bigger seat. The frame is also excellent for use as a mounting frame for molded and other seating systems.
- Uses the same motors and front wheel drive technology as the full size Chairman.
- Programmable through the electronically advanced P&G Pilot+™ controller.
- You can add the Omni+™ alternative drive switch options for a complete high end power chair in a small package.
- Unique elevator/seat tilt feature, allowing active to passive seat positioning, all without tools.

ADDITIONAL FUNCTIONS

- Retractable Joystick
- P&G Pilot+™ and Omni+™ alternative drive switch options
- Wired Remote control panel modification for simple use of a remote joystick
- Tubular interface kit for mounting aftermarket seating systems

OPTIONAL ACCESSORIES

- Foot straps
- 4-point belt
- Seat belt
- Head rest
- Thigh supports
- Lateral supports
- Lexan tray



Retractable joystick



Brake release



Fabric or leatherette



Permobil Inc.

6961 Eastgate Boulevard
Lebanon, TN 37090
Phone 800-736-0925
Fax 800-231-3256

¹ MK battery phone # 1-800-372-9253
² Depending on type of batteries and their condition
Specifications are subject to change.
Does not include custom sizes.
See alternative drive order form for more options.